

ARTICLE

Open Access

Mass spectrometry for the evaluation of monoclonal proteins in multiple myeloma and related disorders: an International Myeloma Working Group Mass Spectrometry Committee Report

David L. Murray¹, Noemi Puig², Sigurdur Kristinsson³, Saad Z. Usmani⁴, Angela Dispenzieri^{1,5}, Giada Bianchi⁶, Shaji Kumar⁵, Wee Joo Chng^{7,8,9}, Roman Hajek¹⁰, Bruno Paiva¹¹, Anders Waage^{12,13}, S. Vincent Rajkumar⁵ and Brian Durie¹⁴

Abstract

Plasma cell disorders (PCDs) are identified in the clinical lab by detecting the monoclonal immunoglobulin (M-protein) which they produce. Traditionally, serum protein electrophoresis methods have been utilized to detect and isotype M-proteins. Increasing demands to detect low-level disease and new therapeutic monoclonal immunoglobulin treatments have stretched the electrophoretic methods to their analytical limits. Newer techniques based on mass spectrometry (MS) are emerging which have improved clinical and analytical performance. MS is gaining traction into clinical laboratories, and has replaced immunofixation electrophoresis (IFE) in routine practice at one institution. The International Myeloma Working Group (IMWG) Mass Spectrometry Committee reviewed the literature in order to summarize current data and to make recommendations regarding the role of mass spectrometric methods in diagnosing and monitoring patients with myeloma and related disorders. Current literature demonstrates that immune-enrichment of immunoglobulins coupled to intact light chain MALDI-TOF MS has clinical characteristics equivalent in performance to IFE with added benefits of detecting additional risk factors for PCDs, differentiating M-protein from therapeutic antibodies, and is a suitable replacement for IFE for diagnosing and monitoring multiple myeloma and related PCDs. In this paper we discuss the IMWG recommendations for the use of MS in PCDs.

Background

Plasma cell disorders (PCDs) are a group of diseases characterized by clonal expansion of plasma cells¹. Central to the diagnosis and monitoring of most PCDs is detection of the monoclonal immunoglobulin

components which are generally overproduced by the expanding plasma cell clone. This overproduced monoclonal immunoglobulin (often referred to as an M-protein or paraprotein) typically is an intact immunoglobulin, and also can be either the free light chain (LC) component alone or the heavy chain component alone in rare instances².

While the M-protein is homogeneous and typically constant in any particular patient, the heterogeneity of M-proteins from patient to patient is significant and thus a diverse set of methods are employed to characterize and

Correspondence: David L. Murray (murray.david@mayo.edu)

¹Department of Laboratory Medicine and Pathology, Mayo Clinic, Rochester, MN, USA

²Hospital Universitario de Salamanca, Instituto de Investigación Biomédica de Salamanca, Salamanca, Spain

Full list of author information is available at the end of the article

© The Author(s) 2021



Open Access This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons license, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons license and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this license, visit <http://creativecommons.org/licenses/by/4.0/>.

detect M-proteins³. Serum protein electrophoresis (SPEP) enables the detection and relative quantitation of the M-protein, whereas serum immunofixation electrophoresis (IFE) enables establishment of M-protein isotype. Another widely utilized assay is the serum free light chain (sFLC) assay that utilizes specific antibodies for quantitation of circulating free kappa (κ) and lambda (λ) light chains (LCs)⁴. The sFLC assay is an analytically sensitive assay for M-protein detection whereby an abnormal κ/λ FLC ratio ($\kappa/\lambda < 0.26$ or > 1.65) suggests the presence of an aberrant plasma or B-cell clone, but not all patients with multiple myeloma (MM) have abnormal sFLC ratios at diagnosis⁵. The sFLC assay has demonstrated a particularly important supportive role in the diagnosis, prognosis, and monitoring of PCDs, especially for amyloidosis and non-secretory MM patients^{1,6,7}. In 2009, a panel of members of the IMWG met to develop guidelines for standard investigative workup of patients with suspected MM. The group evaluated studies that compared the diagnostic sensitivity and specificity of different combinations of testing panels and concluded that a combination of SPEP, serum IFE, and sFLC should be used for screening⁸. These protein studies have been incorporated in defining the various PCDs, as well as the uniform response criteria to assess treatment efficacy⁹.

While the clinical utility of the electrophoretic methods to screen and monitor PCDs has been well-established, several changes in the treatment of PCDs are pushing these methods to their analytical limits. Dramatic improvement in the treatment response of MM patients to new chemo and immune therapies is challenging long-held assumptions about this ominous disease. In particular, as durable remissions can now be achieved for most MM patients, highly sensitive methods capable of detecting low-level disease are important for patient counseling. The long-standing routine serum electrophoretic methods (SPEP and IFE) are not capable of providing the analytical sensitivity needed to assess minimal residual disease (MRD). As a result, laboratories have turned to bone marrow (BM) aspirates and/or biopsies and detecting residual malignant plasma cells with high-sensitivity flow cytometry¹⁰ and their specific VDJH/DJH sequences by next-generation sequencing (NGS)¹¹. These BM-based techniques are well-established and currently available for MRD testing after therapy. However, these MRD methods require an invasive procedure and a lab with a higher level of expertise to perform the testing. Attesting to its impact on prognosis, NGS is now an FDA-approved method for measuring MRD in MM. In addition, new monoclonal therapeutic antibodies (t-mAbs) designed to eradicate malignant plasma cells are producing interferences making it difficult to distinguish between a patient's M-protein and the t-mAb drug on electrophoretic methods¹².

A search for a more convenient serum-based test to complement BM MRD detection and aid in resolving t-mAb interferences was sought to address limitations in traditional testing. Mass spectrometry (MS) is aptly suited for this task as the improvements in instrumentation and techniques have resulted in increased resolution and sensitivity which have outpaced improvements in electrophoresis.

M-protein detection by mass spectrometry

The basis of all the mass spectrometric methods for detecting M-proteins is the unique sequence of the antigen binding region, also called the complementarity determining region (CDR) of the Ig¹³. Each plasma cell produces a unique Ig with a specific CDR due to the adaptive immune system's optimization of the CDR via somatic hypermutation to increase the affinity of the Ig to its target antigen. The resulting CDR amino acid sequence is unique, and each plasma cell clone has a different peptide sequence and overall mass, which is the basis of the M-protein detection by MS¹⁴. Two MS methods have emerged in the literature. Both methods start with immune-enrichment of patient immunoglobulins (Igs) but differ on the analytical target used to detect the M-protein. One method utilizes Ig trypsin digestion and detection of peptides specific to the M-protein CDR^{15–18}. This method has been termed the "clonotypic peptide" approach. The second method utilizes total LC mass distributions from Igs which have been chemically reduced and denatured into heavy and light chain components^{19,20}. This method will be termed intact LC mass measurements herein.

The clonotypic peptide approach

The clonotypic peptide approach has been employed to monitor MM patients¹⁵, assess for MRD in complete response (CR) patients^{16,21}, and to differentiate M-proteins from Daratumumab¹⁷. The clonotypic method is very analytically sensitive with M-protein detection rates down to 0.001 g/L giving it potential to be a serum-based MRD method¹⁶. However, the method has some limitations that hinder it from being suitable for all patients. In some patients with M-proteins above 10 g/L, only constant/framework peptides were detected making it impossible to follow or detect M-proteins in these patients. In addition, as the CDR also contains framework regions, which have less diversity across the polyclonal background, clonotypic peptide selection needs to be done judiciously to assure the uniqueness of the clonotypic peptide to the M-protein. As such, the success of identifying a clonotypic peptide still requires genetic CDR sequence information from the malignant plasma cells. Once a potential clonotypic peptide is determined, a search of the human proteome is necessary to assure its

uniqueness. Unfortunately, a complete library of the possible human CDR tryptic peptides is not yet available. All this requires very thorough understanding of immunoglobulin CDR region to assure definitive test results. Lastly, screening for M-proteins using the clonotypic peptide method has not been demonstrated. Recent work utilizing IFE gels followed by the clonotypic approach has aided in simplifying the pre-analytical work for the clonotypic method²², but does not remove IFE from its screening role.

More clinical data is still needed to determine if the MRD positivity by the clonotypic approach translates into disease-free and overall survival benefits. While most studies utilizing the clonotypic method detect the presence of the M-protein in samples from patients who are BM MRD-negative, the clinical implications of this will require more time to elucidate.

Intact Ig light chain approach

A second more simpler and practical approach for M-protein detection relies on scanning the overall mass distribution of denatured intact Ig LCs. In this method, the polyclonal Ig LC mass distribution results in two normal distributions, one for the lambda and second for the kappa LCs²³. The first version of this intact LC mass measurement method was performed on higher performance mass spectrometers (ESI-Q-TOFs) with liquid chromatography. This method (termed monoclonal immunoglobulin rapid accurate mass measurements or miRAMM) was found to have an analytical limit of quantitation of 0.05 g/L and a limit of detection of 0.01 g/L. Later this method was adapted to matrix-assisted laser desorption/ionization-time-of-flight (MALDI-TOF) mass spectrometry (MS)²⁴, which eliminated the chromatography step and reduced analytical time from 20 min to 10 s. Eventually, the MALDI-TOF method was modified to include immuno-enrichments to conclusively identify the M-protein isotype in a similar fashion to IFE^{20,25,26}. The discussion and recommendations for MS contained in this working group report, unless otherwise specified, pertain to this intact LC MALDI-TOF MS method.

The intact LC MALDI-TOF MS assay is referred to as Mass-Fix at the Mayo Clinic to emphasize its similarity to conventional immunofixation. It has been clinically validated as a lab-developed test at the Mayo Clinic²⁷. At this point in time, over 40,000 clinical samples have been analyzed by the intact LC MALDI-TOF method at Mayo Clinic. The Mass-Fix assay has been found to be comparable to IFE in terms of assay turnaround time and ease of interpretation with the added benefit of reduced labor. These features make the intact LC method more economically attractive in comparison to the clonotypic method. Commercial efforts are being made to automate the technique and provide a high-throughput method

which would be widely accessible for routine laboratory implementation. Validation is currently being performed in several clinical laboratories.

Analytical and clinical performance of intact LC MALDI-TOF assay compared to IFE

Several comparisons of the intact LC MALDI-TOF MS and IFE have been published^{20,21,25}. The ability of intact LC MALDI-TOF MS method to detect M-proteins was retrospectively evaluated on a large cohort of patient samples with physician-ordered SPEP and IFE, which included a wide variety of M-proteins encountered in clinical practice. MALDI-TOF MS detected M-proteins in 100% of samples that were positive by both SPEP and IFE ($n = 84$) and in 97% of samples positive by IFE but negative by SPEP (68 of 70). In samples negative by both SPEP and IFE ($n = 28$), MALDI-TOF identified one positive case²⁰. To test analytical sensitivity, a cohort of 27 PCD patients was diluted into normal human serum and both IFE and intact LC MALDI-TOF were performed. The percentage of M-proteins detected by the intact LC MALDI-TOF was significantly higher than IFE at every dilution. A second prospective study involving 290 patients with documented PCDs demonstrated that substitution of serum and urine intact LC MALDI-TOF testing results for convention IFE results allowed for detection of disease in 18 more patients. Overall, MS had higher clinical sensitivity than IFE in detecting M-proteins in patients with abnormal FLC ratios. A combination of FLC, serum and urine Mass-Fix performed equally as well as current IMWG recommendations of PEL, IFE, and sFLC²⁸. Incorporation of serum FLC MALDI-TOF MS²⁶ has been shown to improve the sensitivity of M-protein detection during monitoring of AL amyloidosis^{29,30}. The intact LC method has also been used for the diagnosis of heavy chain diseases³¹.

A more recent study confirmed the increased sensitivity and accuracy of intact LC MALDI-TOF MS over IFE. In this study, 226 patients from the Olmsted screening cohort³² who were initially negative for a PCD by SPEP but who were found to have a PCD on later follow-up were studied. Since the initial screen was performed by SPEP, the study assessed the ability of IFE, the intact LC MALDI-TOF method (Mass-Fix), as well as a more sensitive miRAMM method to detect the M-protein in the original sample. M-protein would have been detected in 11%, 50%, and 67% of these participants if initial screening was performed by IFE, MALDI-TOF, and miRAMM, respectively³³. A preliminary study of the pre-commercial intact LC MALDI-TOF method (QIP-MS, The Binding Site) utilizing samples from the GEM-CESAT trial patients demonstrated moderate correlation with MRD assessment by next-generation flow MRD status³⁴.

Table 1 List of IMWG recommendations regarding mass spectrometry.

Intact LC MALDI-TOF can be used in lieu of immunofixation in the clinical assessment of patients and the assessment of patients on clinical trials.

We endorse the use of mass spectrometry to aid in distinguishing therapeutic antibodies from endogenous M-proteins.

We recognize that using mass spectrometric methods in lieu of conventional IFE may lead to lower rates of CR, and therefore cross-comparisons of CR rates in trials done in different time periods is not recommended.

We endorse the collection of further data from mass spectrometry (MALDI-TOF, miRAMM, or clonotypic peptide approach) to document the ability to test for MRD negativity in the peripheral blood, and to guide timing of BM tests for next-generation flow cytometry and NGS studies.

We encourage further investigations to further clarify the relationship and implications of N-linked glycosylation in MGUS progression to myeloma and AL amyloidosis.

The intact LC method MS can give useful information to discriminate between t-mAbs and M-proteins as the mass of the t-mAb LC can be measured. Using the intact LC MALDI-TOF method, daratumumab and IgG kappa M-proteins were distinguishable in 84% of the samples tested³⁵. The inability to discriminate between daratumumab and M-protein was most likely due to the lower resolution of the MALDI-TOF assay, as studies involving the higher resolution LC-ESI-QTOF (miRAMM) MS had a 100% success rate of differentiating the t-mAb (daratumumab, elotuzumab, and isatuximab) from M-protein³⁶.

A previously underappreciated finding of M-proteins when employing the intact LC method was the observation of M-proteins with a LC mass that was broad and outside the normal mass ranges of kappa and lambda LCs^{25,28}. Additional studies demonstrated that these LCs contained N-linked glycosylation and were to be more prevalent in patients with AL-amyloidosis³⁷ and cold agglutinin disease³⁸ than other PCDs. Interestingly, LC glycosylation was also found to be an independent risk factor for the progression of MGUS to myeloma and related disorders³⁹. These glycosylated LCs were demonstrated to be present years before the diagnosis of myeloma of related disorder and hence may allow for identification of patients at higher risk for PCDs⁴⁰.

MS using the more sensitive miRAMM intact LC approach may have a role in MRD detection. In one study, residual monoclonal proteins were detected in patients who were in stringent CR (sCR) (10^{-4} to 10^{-5} by 6 color flow cytometry) in 80% ($n = 16$) and 60% ($n = 25$) of patients at 6–12 months, respectively⁴¹. The influence of Ig half-life on the MRD data needs to be further elucidated. However, since miRAMM MS is a serum-based assay, the test could be performed serially allowing the rate of response to be assessed over time. A second more recent study also compared the intact LC MALDI-TOF MS with BM flow cytometry-based MRD⁴². This exploratory study demonstrated MALDI-TOF MS was concordant with MRD flow for 44/71 (62%) patients ($p = 0.342$). Amongst the 27 discordant cases, 17 were

detectable only by MALDI-TOF MS, and 10 detectable only by MRD flow. Additional prospective studies are required to delineate the role and timing of these modalities in assessing MRD.

Recommendations

We conclude that MS has the advantage of increased accuracy, documented clinical and analytic sensitivity, and the intact LC MALDI-TOF method is easier on laboratory work flow for the detection of M-proteins²⁰. As summarized in Table 1, the IMWG Mass Spectrometry Committee endorses detection of M-proteins by MS (intact MALDI-TOF method) as an alternative to IFE for clinical practice and clinical trials. The group also endorses MS for distinguishing residual M-protein from therapeutic monoclonal antibodies for clinical practice, and for accurate interpretation and determination of complete response in clinical trials. We recognize that using mass spectrometric methods instead of conventional IFE may lead to lower rates of CR, and therefore cross-comparisons of CR rates in trials done in different time periods is not recommended. We hope, with further data, that mass spectrometric methods (MALDI-TOF, miRAMM, or clonotypic peptide approach) may provide the ability to test for measurable disease in the peripheral blood and help guide timing of BM tests for next-generation flow cytometry and NGS studies. Prospective validation of the impact of MRD negativity based on MS on patient progression-free survival and overall survival in relationship with BM NGS MRD is strongly encouraged as part of clinical trials. Finally, we encourage further investigations to further clarify the relationship and implications of N-linked glycosylation and AL amyloidosis.

Author details

¹Department of Laboratory Medicine and Pathology, Mayo Clinic, Rochester, MN, USA. ²Hospital Universitario de Salamanca, Instituto de Investigación Biomédica de Salamanca, Salamanca, Spain. ³Faculty of Medicine, University of Iceland, Reykjavik, Iceland. ⁴Department of Hematologic Oncology & Blood Disorders, Levine Cancer Institute/Atrium Health, Charlotte, NC, USA. ⁵Department of Hematology, Mayo Clinic, Rochester, MN, USA. ⁶Division of Hematology, Department of Medicine, Brigham and Women's Hospital,

Harvard Medical School, Boston, MA, USA. ⁷Cancer Science Institute of Singapore, NUS, Singapore, Singapore. ⁸Yong Loo Lin School of Medicine, NUS, Singapore, Singapore. ⁹National University Cancer Institute, Singapore, Singapore. ¹⁰Department of Hematooncology, University Hospital Ostrava and Faculty of Medicine, University of Ostrava, Ostrava, Czech Republic. ¹¹Clinica Universidad de Navarra, Centro de Investigacion Medica Aplicada (CIMA), IDISNA, Pamplona, Spain. ¹²Department of Clinical and Molecular Medicine, Norwegian University of Science and Technology, Trondheim, Norway. ¹³Department of Hematology, St. Olav's University Hospital, Trondheim, Norway. ¹⁴Department of Hematology, Cedars-Sinai Outpatient Cancer Center, Los Angeles, CA, USA

Conflict of interest

D. Murray has intellectual property right with potential royalties that were licensed to The Binding Site.

Publisher's note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Supplementary information The online version contains supplementary material available at <https://doi.org/10.1038/s41408-021-00408-4>.

Received: 24 August 2020 Revised: 10 November 2020 Accepted: 1 December 2020

Published online: 01 February 2021

References

- Rajkumar, S. V. et al. International Myeloma Working Group updated criteria for the diagnosis of multiple myeloma. *Lancet Oncol.* **15**, e538–e548 (2014).
- Kyle, R. A. et al. Review of 1027 patients with newly diagnosed multiple myeloma. *Mayo Clinic Proc.* **78**, 21–33 (2003).
- Katzmann, J. A. et al. Screening panels for detection of monoclonal gammopathies. *Clin. Chem.* **55**, 1517–1522 (2009).
- Bradwell, A. R. et al. Highly sensitive, automated immunoassay for immunoglobulin free light chains in serum and urine. *Clin. Chem.* **47**, 673–680 (2001).
- Katzmann, J. A. et al. Serum reference intervals and diagnostic ranges for free kappa and free lambda immunoglobulin light chains: relative sensitivity for detection of monoclonal light chains. *Clin. Chem.* **48**, 1437–1444 (2002).
- Dispenzieri, A. et al. International Myeloma Working Group guidelines for serum-free light chain analysis in multiple myeloma and related disorders. *Leukemia.* **23**, 215–224 (2009).
- Dejoie, T. et al. Serum free light chains, not urine specimens, should be used to evaluate response in light-chain multiple myeloma. *Blood.* **128**, 2941–2948 (2016).
- Dimopoulos, M. et al. Consensus recommendations for standard investigative workup: report of the International Myeloma Workshop Consensus Panel 3. *Blood.* **117**, 4701–4705 (2011).
- Kumar, S. et al. International Myeloma Working Group consensus criteria for response and minimal residual disease assessment in multiple myeloma. *Lancet Oncol.* **17**, e328–e346 (2016).
- Rawstron, A. C. et al. Minimal residual disease assessed by multiparameter flow cytometry in multiple myeloma: impact on outcome in the Medical Research Council Myeloma IX Study. *J. Clin. Oncol.* **31**, 2540–2547 (2013).
- Martinez-Lopez, J. et al. Prognostic value of deep sequencing method for minimal residual disease detection in multiple myeloma. *Blood.* **123**, 3073–3079 (2014).
- Mills, J. R. & Murray, D. L. Identification of friend or foe: the laboratory challenge of differentiating M-proteins from monoclonal antibody therapies. *J. Appl. Lab. Med.* **1**, 421–431 (2017).
- Zajec, M. et al. Mass spectrometry for identification, monitoring, and minimal residual disease detection of M-proteins. *Clin. Chem.* **66**, 421–433 (2020).
- Dekker, L. J. M. et al. An antibody-based biomarker discovery method by mass spectrometry sequencing of complementarity determining regions. *Anal. Bioanal. Chem.* **399**, 1081–1091 (2010).
- Barnidge, D.R. et al. Monitoring M-proteins in patients with multiple myeloma using heavy-chain variable region clonotypic peptides and LC-MS/MS. *J. Proteome Res.* **13**, 1905–1910 (2014).
- Bergen, H. R. et al. Clonotypic light chain peptides identified for monitoring minimal residual disease in multiple myeloma without bone marrow aspiration. *Clin. Chem.* **62**, 243–251 (2016).
- Zajec, M. et al. Development of a targeted mass-spectrometry serum assay to quantify M-protein in the presence of therapeutic monoclonal antibodies. *J. Proteome Res.* **17**, 1326–1333 (2018).
- Remily-Wood, E. R. et al. Quantification of peptides from immunoglobulin constant and variable regions by LC-MRM MS for assessment of multiple myeloma patients. *Proteomics Clin. Appl.* **8**, 783–795 (2014).
- Barnidge, D.R. et al. Using mass spectrometry to monitor monoclonal immunoglobulins in patients with a monoclonal gammopathy. *J. Proteome Res.* **13**, 1419–1427 (2014).
- Mills, J. R. et al. Comprehensive assessment of M-proteins using nanobody enrichment coupled to MALDI-TOF mass spectrometry. *Clin. Chem.* **62**, 1334–1344 (2016).
- Martins, C. O. et al. Mass spectrometry-based method targeting Ig variable regions for assessment of minimal residual disease in multiple myeloma. *J. Mol. Diagn.* **22**, 901–911 (2020).
- Zajec, M. et al. Integrating serum protein electrophoresis with mass spectrometry, a new workflow for M-protein detection and quantification. *J. Proteome Res.* **19**, 2845–2853 (2020).
- Barnidge, D. R. et al. Phenotyping polyclonal kappa and lambda light chain molecular mass distributions in patient serum using mass spectrometry. *J. Proteome Res.* **13**, 5198–5205 (2014).
- Barnidge, D. R., Griffin, T. J. & Murray, D. L. Using matrix-assisted laser desorption/ionization time-of-flight mass spectrometry to detect monoclonal immunoglobulin light chains in serum and urine. *Rapid Commun. Mass Spectrom.* **29**, 1–4 (2015).
- Kohlhagen, M. C. et al. Screening method for M-proteins in serum using nanobody enrichment coupled to MALDI-TOF mass spectrometry. *Clin. Chem.* **62**, 1345–1352 (2016).
- Sepiashvili, L. et al. Direct detection of monoclonal free light chains in serum by use of immunoenrichment-coupled MALDI-TOF mass spectrometry. *Clin. Chem.* **65**, 1015–1022 (2019).
- Kohlhagen, M.C. et al. Automation and validation of a MALDI-TOF MS (Mass-Fix) replacement of immunofixation electrophoresis in the clinical lab. *Clin. Chem. Lab. Med.* **59**, 155–163 (2020).
- Milani, P. et al. The utility of MASS-FIX to detect and monitor monoclonal proteins in the clinic. *Am. J. Hematol.* **92**, 772–779 (2017).
- Sharpley, F. A. et al. A novel mass spectrometry method to identify the serum monoclonal light chain component in systemic light chain amyloidosis. *Blood Cancer J.* **9**, 16 (2019).
- Dispenzieri, A. et al. Blood mass spectrometry detects residual disease better than standard techniques in light-chain amyloidosis. *Blood Cancer J.* **10**, 20 (2020).
- Thoren, K. L. et al. Identification of gamma heavy chain disease using MALDI-TOF mass spectrometry. *Clin. Biochem.* **77**, 57–61 (2020).
- Kyle, R. A. et al. A long-term study of prognosis in monoclonal gammopathy of undetermined significance. *N. Engl. J. Med.* **346**, 564–569 (2002).
- Murray, D. L. et al. Detection and prevalence of monoclonal gammopathy of undetermined significance: a study utilizing mass spectrometry-based monoclonal immunoglobulin rapid accurate mass measurement. *Blood Cancer J.* **9**, 102 (2019).
- Puig, N. et al. Qip-mass spectrometry in high risk smoldering multiple myeloma patients included in the GEM-CESAR Trial: comparison with conventional and minimal residual disease IMWG response assessment. *Blood* **134**, 581 (2019).
- Moore, L. M., Cho, S. & Thoren, K. L. MALDI-TOF mass spectrometry distinguishes daratumumab from M-proteins. *Clin. Chim. Acta* **492**, 91–94 (2019).
- Mills, J. R. et al. A universal solution for eliminating false positives in myeloma due to therapeutic monoclonal antibody interference. *Blood* **132**, 670–672 (2018).

37. Kumar, S. et al. Assay to rapidly screen for immunoglobulin light chain glycosylation: a potential path to earlier AL diagnosis for a subset of patients. *Leukemia* **33**, 254–257 (2019).
38. Sidana, S. et al. Glycosylation of immunoglobulin light chains is highly prevalent in cold agglutinin disease. *Am. J. Hematol.* **95**, E222–E225 (2020).
39. Dispenzieri, A. et al. N-glycosylation of monoclonal light chains on routine MASS-FIX testing is a risk factor for MGUS progression. *Leukemia* **34**, 2749–2753 (2020).
40. Kourelis, T. et al. MASS-FIX may allow identification of patients at risk for light chain amyloidosis before the onset of symptoms. *Am. J. Hematol.* **93**, E368–E370 (2018).
41. Mills, J. R., Barnidge, D. R., Dispenzieri, A. & Murray, D. L. High sensitivity blood-based M-protein detection in sCR patients with multiple myeloma. *Blood Cancer J.* **7**, e590 (2017).
42. Eveillard, M. et al. Comparison of MALDI-TOF mass spectrometry analysis of peripheral blood and bone marrow-based flow cytometry for tracking measurable residual disease in patients with multiple myeloma. *Br. J. Haematol.* **189**, 904–907 (2020).

Terms and Conditions

Springer Nature journal content, brought to you courtesy of Springer Nature Customer Service Center GmbH (“Springer Nature”).

Springer Nature supports a reasonable amount of sharing of research papers by authors, subscribers and authorised users (“Users”), for small-scale personal, non-commercial use provided that all copyright, trade and service marks and other proprietary notices are maintained. By accessing, sharing, receiving or otherwise using the Springer Nature journal content you agree to these terms of use (“Terms”). For these purposes, Springer Nature considers academic use (by researchers and students) to be non-commercial.

These Terms are supplementary and will apply in addition to any applicable website terms and conditions, a relevant site licence or a personal subscription. These Terms will prevail over any conflict or ambiguity with regards to the relevant terms, a site licence or a personal subscription (to the extent of the conflict or ambiguity only). For Creative Commons-licensed articles, the terms of the Creative Commons license used will apply.

We collect and use personal data to provide access to the Springer Nature journal content. We may also use these personal data internally within ResearchGate and Springer Nature and as agreed share it, in an anonymised way, for purposes of tracking, analysis and reporting. We will not otherwise disclose your personal data outside the ResearchGate or the Springer Nature group of companies unless we have your permission as detailed in the Privacy Policy.

While Users may use the Springer Nature journal content for small scale, personal non-commercial use, it is important to note that Users may not:

1. use such content for the purpose of providing other users with access on a regular or large scale basis or as a means to circumvent access control;
2. use such content where to do so would be considered a criminal or statutory offence in any jurisdiction, or gives rise to civil liability, or is otherwise unlawful;
3. falsely or misleadingly imply or suggest endorsement, approval, sponsorship, or association unless explicitly agreed to by Springer Nature in writing;
4. use bots or other automated methods to access the content or redirect messages
5. override any security feature or exclusionary protocol; or
6. share the content in order to create substitute for Springer Nature products or services or a systematic database of Springer Nature journal content.

In line with the restriction against commercial use, Springer Nature does not permit the creation of a product or service that creates revenue, royalties, rent or income from our content or its inclusion as part of a paid for service or for other commercial gain. Springer Nature journal content cannot be used for inter-library loans and librarians may not upload Springer Nature journal content on a large scale into their, or any other, institutional repository.

These terms of use are reviewed regularly and may be amended at any time. Springer Nature is not obligated to publish any information or content on this website and may remove it or features or functionality at our sole discretion, at any time with or without notice. Springer Nature may revoke this licence to you at any time and remove access to any copies of the Springer Nature journal content which have been saved.

To the fullest extent permitted by law, Springer Nature makes no warranties, representations or guarantees to Users, either express or implied with respect to the Springer nature journal content and all parties disclaim and waive any implied warranties or warranties imposed by law, including merchantability or fitness for any particular purpose.

Please note that these rights do not automatically extend to content, data or other material published by Springer Nature that may be licensed from third parties.

If you would like to use or distribute our Springer Nature journal content to a wider audience or on a regular basis or in any other manner not expressly permitted by these Terms, please contact Springer Nature at

onlineservice@springernature.com